



The following *Spectris* MR Injector component names and acronyms are used interchangeably in this manual:

1. HSU (Injector Head/Stand Unit)
2. HCU (Head Control Unit)
3. OTR (Optical Transmitter/Receiver)
4. DCU (Display Control Unit)
5. MCU (Master Control Unit) in tower configuration
6. HSU (Head/Stand Unit) in mobile configuration

NOTE: Position the HCU in the scanner room so that the HCU is placed in a location where the magnetic field does not exceed 800 gauss within 1 meter of the magnet. To ensure proper placement and operation of the HCU, follow the guidelines in the *Spectris* MR Installation Manual and refer to the site copy of the magnetic fringe field plot provided by the scanner manufacturer. The HCU must be located on the floor.

NOTE: Follow all institutional, local or national safety regulations related to running cables across the floor.

SanMedical GmbH

Gewerbegebiet

Brand No. 55

D-91729 Haundorf

Germany HRB:4520

fon: +49(0)173 576 25 76

email: sanmedical@email.de